



Hop Latent Viroid Levels and Distribution in Cannabis Plant Tissue

Jama, Aisha; Kaprikorn Nursery; Garrido-Lecca, Alfonso; Saldi, Tassa

Table of contents

| | |
|------------------------|----|
| Summary | |
| Overview | 2 |
| Results and Discussion | |
| <hr/> | |
| Follow-up Study | 5 |
| <hr/> | |
| Recommendations | |
| Limitations | 6 |
| References | |
| <hr/> | |
| Supplementary Data | 7 |
| <hr/> | |
| Contact the TUMI Team | 12 |
| <hr/> | |

Summary

The purpose of this study was to determine the variability of Hop Latent Viroid (HLVd) load throughout infected cannabis plants as well as the tissue-type that most reliably predicts infection. While the number of plants assayed in this preliminary study was limited, the data presented below strongly suggest that upper root tissue has the highest viroid load and shows consistent detection of HLVd genetic material in infected cannabis plants.

Overview

Four cannabis plants (teenage) identified as infected with HLVd were acquired by Kaprikorn Nurseries for the purposes of this study. Infected plants were kept in strict isolation. To determine variability in viroid amounts throughout the plants, approximately twenty-five individual samples were taken from distinct locations around each plant. Tissue types sampled include: petiole samples from older growth, petiole samples from newer growth, stem sample from freshly cut upper clones, upper root mass samples and lower root mass samples. A total of 95 samples were assayed in this study. A subset of samples taken from two non-infected plants were also included in this study to confirm the absence of cross-contamination.

All samples were taken within the same three-day time span and were placed in identical TUMI Genomics collection tubes. Samples were mailed and processed in parallel at the TUMI Genomics Laboratory using a triplex TaqMan qRT-PCR assay with a limit of detection (i.e. the viroid level at which the PCR test is 95% accurate) of 7.5 viroid copies per reaction. Additional information about this assay and a full validation report can be found here: [HLVd_validation](#)). All samples were assayed for an internal plant RNA control to confirm sample integrity and alongside three external controls: a non-infected control, an infected control, and a nontemplate control; to rule out assay dysfunction or contamination.

Results and Discussion

Root tissue gives the most consistent HLVd amplification in individual plants. Among the individual samples taken from the four test plants, upper foliage samples showed variable success for detection of HLVd (Table 1). Petiole taken from new growth showed between 10% to 100% success depending on the plant with an average success rate of 36% (9/25). Old growth and tissue from freshly cut clones showed slightly better detection success with an average of 75% (19/25) and 61% (11/18) testing positive for HLVd, respectively. Interestingly, root tissue showed an impressive level of accuracy.

Every root tissue sample taken from either the upper or lower portion of the root ball showed robust amplification of HLVd (36/36 samples testing positive). We saw no HLVd amplification from samples taken from plants with no known HLVd infection (Table 2). These results indicate that inclusion of root tissue when testing for a potential HLVd infection could improve detection accuracy.

We did note that 20 of the samples testing positive in the upper foliage were defined as lowlevel positive (viroid load <15 HLVD copies per reaction). This observation was prominent in both Plant 2 and Plant 3 where almost all positive samples taken from upper foliage showed very low viroid load, suggesting that detection of HLVD in foliage tissue from these plants would likely be missed by lower sensitivity assays. While it is reported that viroids reach extremely high copy number during infection, our results suggest that not all tissues in infected cannabis plants contain high levels of HLVD. Tissues with low viroid load may contain viroid that is trafficking through the vascular system or viroid that has just gained access to the cellular environment, but is not yet efficiently replicating.

Once a system-wide infection is established, the location where the sample is taken from is likely less critical. In support of this, one of the tested plants (Plant 4) showed consistent, robust HLVD amplification in all tissues sampled suggesting development of a full systemic infection in this plant.





Figure 1. (above) Root tissue contains consistent levels of HLVD in positive plants. The percent of samples from each of the four assayed plants where HLVD was detected in the indicated tissue. For each tissue-type, five individual samples were collected from each plant.

Figure 2. (right) The percent success of HLVD detection in indicated tissues across all plants tested. The red box highlights the root tissue, where the detection level of HLVD genetic material was 100%.



Upper root tissue shows higher, less variable levels of HLVd. Because commercial assays detect HLVd with variable sensitivity, understanding the level of HLVd accumulation in different tissues is helpful to estimate the success rate of HLVd detection across different tests. Using standard curve data generated from known concentrations of HLVd sequence (done in triplicate) we estimated the median viroid load across positive samples from each tissue type (Figure 2).

Due to the low number of new growth samples with detectable levels of HLVd, we did not include this tissue in the analysis. As shown in Figure 2, we found that upper root tissue had substantially higher viroid load than other plant tissue analyzed. Upper root tissue had an average of ~500-fold to ~2000-fold more HLVd than foliage tissue analyzed in this study. Even when only considering samples taken from the root, upper root tissue amplified an average of 36-fold more HLVd than samples taken from the lower root.

We also noted that root tissue (both upper and lower) showed less variability in viroid load both within each individual plant and when all tissue samples of each type were considered together. These results indicate that upper root tissue accumulates high-levels of viroid titer and therefore inclusion of upper root tissue when screening for HLVd infection could improve detection sensitivity.





Figure 3. (above) Upper Root tissues accumulates high levels of HLVd. A. Median viroid load and variation across the indicated tissue types. The dashed red line shows the threshold of reliable (>95%) HLVd detection for the TUMI Genomics HLVd PCR assay. New growth samples were not included in this analysis due to a small number of positive samples obtained from this tissue.

Figure 4. (right) Model of viroid trafficking from the site of inoculation to “sink” plant tissues. Image shows two potential sites of viroid introduction (red and blue circles). Note that the location the viroid travels in the foliage is dependent on where the viroid was introduced in the plant. However, HLVd will always travel down to the root system regardless of where HLVd was introduced. Image adapted from Takeda and Ding, 2009.



These findings are consistent with previous studies showing that viroids accumulate to high levels in root tissue, including HLVd in Hops plants (Matousek et. al. 1995, Palukaitis, 1987, Lin et. al, 2015, Nabeshima et. al, 2017, Antignus, et. al. 2007). Because viruses and viroids traffic through the vascular system of plants, they traffic from the site of introduction to areas of active growth, i.e. growth sinks. Because the location of viroid introduction is not usually known, it is difficult to determine which leaves and petioles represent growth sinks at the time of initial infection and are therefore likely to contain spreading viroid. However, the root system is always growing and receiving energy from the upper portion of the plant, so viroids from the site of inoculation will consistently traffic to the root system (Figure 2A) (Takeda and Ding, 2009).

Follow-up Study

To confirm the findings described above, we performed a follow-up study using additional cannabis varieties collected from a distinct, off-site location. Individual samples from sixteen plants of unknown HLVd infection status we collected. For each plant, two types of samples were prepared: 1) three petiole samples from foliage along with a single pooled sample from the same petiole tissues and 2) three root samples along with a single pooled root sample from the same root tissue. Samples collected were within the same day, placed in collection tubes and mailed to TUMI Genomics where they were tested for the presence of HLVd using our RTPCR assay.

Among the tested plants, three out of the sixteen tested positive for HLVd. The success rate of HLVd detection in each root sample and the pooled root samples remained 100%, replicating the data presented above. The foliage samples again showed poor individual success at detecting HLVd infection (Table 1). On one plant, HLVd was not detected in any of the petiole sample but was still present in all three of the individual root samples and the pooled samples.

Table 1. (below) Success rate of HLVd detection among samples taken from positive plants

| Sample | Individual Petiole Success Rate | Pooled Petiole Result | Individual Root Success Rate | Pooled Root Result |
|---------|---------------------------------|-----------------------|------------------------------|--------------------|
| Plant 1 | 33% | Positive | 100% | Positive |
| Plant 2 | 0% | Negative | 100% | Positive |
| Plant 3 | 66% | Positive | 100% | Positive |

We also noted that the TUMI Genomics test was able to detect HLVd in pooled samples without a loss in sensitivity (see supplemental tables), allowing inclusion of multiple samples from a single root mass which can increase sensitivity and accuracy of the test results.

Recommendations

HLVd can exist at dramatically different levels across different tissues within an infected plant and data shown here confirms this assertion. For reliable diagnostic testing, it is best to collect multiple tissue samples from around the plant. This ensures that if the infection is early and/or local (not spread system-wide), diagnostic assays can still detect HLVd in some of the submitted tissue. Based on the data presented here, along with findings in published literature from different viroids in other crops, we recommend that upper root tissue be included when submitting tissue samples for HLVd testing. Because of the limited sample size of this study, we recommend growers continue to also include petiole samples from older growth among the tissue submitted for testing. If root tissue is not available (such as in a young clone), care should be taken to submit multiple tissue samples from older growth areas of the plant.

Limitations

While this study strongly suggests that root tissue is the most reliable tissue to identify HLVd infection, the sample size tested was limited. Seven plants were assayed originating from five different varieties. We cannot completely rule out the possibility that the viroid distribution in the tested plants is cultivar specific. However, given the 100% detection rate of hop latent viroid in every root sample assayed (45 individual samples total), the data clearly indicates that inclusion of root tissue when testing for HLVd can dramatically improve identification of true positive plants.

References

Antignus, Y., Lachman, O., and Pearlsman, M. The spread of Tomato apical stunt viroid (TASVd) in greenhouse tomato crops is associated with seed transmission and bumble bee activity. *Plant Dis.*, 2007 91:47-50

Chun-Yi Lin, Meng-Ling Wu, Tang-Long Shen & Ting-Hsuan Hung. A mutual titer-enhancing relationship and similar localization patterns between Citrus exocortis viroid and Hop stunt viroid co-infecting two citrus cultivars. *Virology Journal* volume 12, 2015, Article number: 142

J Matousek , L Trněná, P Svoboda, P Oriniaková, C P Lichtenstein. The gradual reduction of viroid levels in hop mericlones following heat therapy: a possible role for a nuclease degrading dsRNA. *Biol Chem Hoppe Seyler*, 1995 Dec;376(12):715-21

Palukaitis, Peter. Potato Spindle Tuber Viroid: Investigation of the long-distance, intra-plant transport route. *Virology*, 1987 153, 239-241

Ryuta Takeda and Biao Ding. Viroid Intercellular Trafficking: RNA Motifs, Cellular Factors and Broad Impacts. *Viruses* 2009, 1, 210-221

Tomoyuki Nabeshima, Motoaki Doi, and Munetaka Hosokawa. Comparative Analysis of Chrysanthemum Stunt Viroid Accumulation and Movement in Two Chrysanthemum (*Chrysanthemum morifolium*) Cultivars with Differential Susceptibility to the Viroid Infection. *Front Plant Sci.* 2017; 8: 1940

Supplementary Data

Tables Show the cycle threshold values and estimated viroid load obtained for all samples tested in this study. Samples that tested at average to high positive levels are indicated in red. Samples that tested low level positive (<15 viroid copies per microliter) are shown in blue. NaN indicates non-detectable. Run controls are shown at the end of Table 2.

Table 2. (below) Amplification results among samples taken from HLVd-infected plants

| Plant | Tissue | Result | Control CT Value | HLVd Target 1 CT Value | HLVd Target 2 CT Value | Viroid Load (copies/μl) |
|---------|-------------|--------------|------------------|------------------------|------------------------|-------------------------|
| Plant 1 | New Petiole | Positive | 30.16 | 18.92 | 21.69 | 88090.82 |
| Plant 1 | New Petiole | Positive | 34 | 18.28 | 20.86 | 138731.51 |
| Plant 1 | New Petiole | Negative | 33.77 | NaN | NaN | Not Detected |
| Plant 1 | New Petiole | Negative | 33.97 | NaN | NaN | Not Detected |
| Plant 1 | New Petiole | Negative | 34.16 | NaN | NaN | Not Detected |
| Plant 2 | New Petiole | Low Positive | 31.81 | 40.6 | NaN | 0.02 |
| Plant 2 | New Petiole | Negative | 32.76 | NaN | NaN | Not Detected |
| Plant 2 | New Petiole | Negative | 32.34 | NaN | NaN | Not Detected |
| Plant 2 | New Petiole | Negative | 31.84 | NaN | NaN | Not Detected |
| Plant 2 | New Petiole | Negative | 31.23 | NaN | NaN | Not Detected |
| Plant 3 | New Petiole | Low Positive | 31.14 | 39.03 | NaN | 0.06 |
| Plant 3 | New Petiole | Negative | 32.57 | NaN | NaN | Not Detected |
| Plant 3 | New Petiole | Negative | 30.63 | NaN | NaN | Not Detected |
| Plant 3 | New Petiole | Negative | 30.22 | NaN | NaN | Not Detected |
| Plant 3 | New Petiole | Negative | 30.4 | NaN | NaN | Not Detected |
| Plant 4 | New Petiole | Positive | 30.36 | 18.05 | 21.05 | 163327.7 |
| Plant 4 | New Petiole | Positive | 30.38 | 19.63 | 22.59 | 53224.62 |
| Plant 4 | New Petiole | Positive | 30.45 | 21.45 | 24.59 | 14628.46 |
| Plant 4 | New Petiole | Positive | 30.18 | 18.73 | 21.7 | 100806.28 |
| Plant 4 | New Petiole | Low Positive | 30.25 | 35.7 | NaN | 0.59 |
| Plant 1 | Old Petiole | Positive | 30.1 | 24.19 | 26.83 | 2092.89 |
| Plant 1 | Old Petiole | Positive | 33.07 | 27.42 | 30.44 | 211.48 |
| Plant 1 | Old Petiole | Positive | 30.21 | 20.52 | 24.15 | 28302.15 |
| Plant 1 | Old Petiole | Positive | 29.35 | 26.68 | 27.16 | 1478.19 |
| Plant 1 | Old Petiole | Negative | 29.44 | NaN | NaN | Not Detected |
| Plant 2 | Old Petiole | Negative | 28.17 | NaN | NaN | Not Detected |
| Plant 2 | Old Petiole | Positive | 28.36 | 29.27 | 35.41 | 56.90 |
| Plant 2 | Old Petiole | Low Positive | 28.17 | 38.03 | NaN | 0.11 |
| Plant 2 | Old Petiole | Low Positive | 28.21 | 36.71 | NaN | 0.29 |
| Plant 2 | Old Petiole | Low Positive | 28.63 | 33.32 | NaN | 3.21 |

| Plant | Tissue | Result | Control CT Value | HLVd Target 1 CT Value | HLVd Target 2 CT Value | Viroid Load (copies/μl) |
|---------|-------------|--------------|------------------|------------------------|------------------------|-------------------------|
| Plant 3 | Old Petiole | Low Positive | 28.11 | 31.36 | 36 | 12.91 |
| Plant 3 | Old Petiole | Positive | 28.66 | 28.34 | 33.9 | 110.09 |
| Plant 3 | Old Petiole | Low Positive | 29.13 | 37.7 | NaN | 0.14 |
| Plant 3 | Old Petiole | Low Positive | 29.52 | 35.97 | NaN | 0.49 |
| Plant 3 | Old Petiole | Negative | 29.34 | NaN | NaN | Not Detected |
| Plant 4 | Old Petiole | Positive | 27.94 | 20.63 | 22.55 | 26176.89 |
| Plant 4 | Old Petiole | Positive | 29.6 | 25.52 | 28.52 | 814.42 |
| Plant 4 | Old Petiole | Positive | 28.35 | 17.71 | 21.06 | 207895.93 |
| Plant 4 | Old Petiole | Positive | 29.47 | 21.45 | 24.52 | 15628.46 |
| Plant 4 | Old Petiole | Low Positive | 27.64 | 31.84 | NaN | 9.18 |
| Plant 1 | Lower Root | Positive | 28.49 | 26.15 | 29.76 | 520.82 |
| Plant 1 | Lower Root | Positive | 31.65 | 27.9 | 31.86 | 150.43 |
| Plant 1 | Lower Root | Positive | 27.35 | 21.88 | 25.42 | 10781.39 |
| Plant 1 | Lower Root | Positive | 25.45 | 19.63 | 22.37 | 53224.62 |
| Plant 2 | Lower Root | Positive | 25.75 | 19.77 | 22.45 | 48190.93 |
| Plant 2 | Lower Root | Positive | 24.09 | 18.62 | 21.41 | 108990.59 |
| Plant 2 | Lower Root | Positive | 27 | 23.51 | 26.07 | 3390.92 |
| Plant 2 | Lower Root | Positive | 27.12 | 25.72 | 29.1 | 706.66 |
| Plant 3 | Lower Root | Positive | 26.04 | 26.23 | 28.11 | 492.07 |
| Plant 3 | Lower Root | Positive | 30.27 | 26.98 | 30.14 | 288.90 |
| Plant 3 | Lower Root | Positive | 28.25 | 24.09 | 26.98 | 2246.80 |
| Plant 3 | Lower Root | Positive | 24.7 | 20.4 | 23.54 | 30817.88 |
| Plant 4 | Lower Root | Positive | 28.05 | 27.09 | 29.72 | 267.29 |
| Plant 4 | Lower Root | Positive | 28.26 | 28.06 | 33.5 | 134.29 |
| Plant 4 | Lower Root | Positive | 29.39 | 28.34 | 32.87 | 110.09 |
| Plant 4 | Lower Root | Positive | 31.49 | 27.55 | 30.35 | 192.85 |
| Plant 1 | Upper Root | Positive | 32.48 | 25.81 | 28.33 | 662.93 |
| Plant 1 | Upper Root | Positive | 27.77 | 18.64 | 21.24 | 107454.62 |
| Plant 1 | Upper Root | Positive | 27.31 | 21.15 | 24.45 | 18099.10 |
| Plant 1 | Upper Root | Positive | NaN | 27.58 | 33.22 | 188.78 |
| Plant 1 | Upper Root | Positive | 25.09 | 18.41 | 21.67 | 126505.62 |
| Plant 2 | Upper Root | Positive | 24.48 | 19.36 | 21.95 | 64465.16 |
| Plant 2 | Upper Root | Positive | 23.91 | 19.05 | 22.11 | 80327.70 |
| Plant 2 | Upper Root | Positive | 25.66 | 20.13 | 23.28 | 37326.33 |
| Plant 2 | Upper Root | Positive | 24.48 | 23.64 | 25.58 | 3092.09 |
| Plant 2 | Upper Root | Positive | 26.23 | 20.41 | 23.54 | 30599.96 |

| Plant | Tissue | Result | Control CT Value | HLVd Target 1 CT Value | HLVd Target 2 CT Value | Viroid Load (copies/μl) |
|---------|------------------|--------------|------------------|------------------------|------------------------|-------------------------|
| Plant 3 | Upper Root | Positive | 28.32 | 23.03 | 26.13 | 4767.07 |
| Plant 3 | Upper Root | Positive | 26.09 | 23.02 | 26.07 | 4801.02 |
| Plant 3 | Upper Root | Positive | 24.51 | 20.62 | 22.93 | 26363.31 |
| Plant 3 | Upper Root | Positive | 27.28 | 22.84 | 25.14 | 5455.17 |
| Plant 3 | Upper Root | Positive | 27.54 | 21.02 | 24.3 | 19848.25 |
| Plant 4 | Upper Root | Positive | 27.23 | 23.61 | 26.56 | 3158.63 |
| Plant 4 | Upper Root | Positive | 25.43 | 19.21 | 22.88 | 71705.79 |
| Plant 4 | Upper Root | Positive | 25.28 | 20.24 | 22.81 | 34523.43 |
| Plant 4 | Upper Root | Positive | 28.65 | 25.12 | 27.69 | 1081.74 |
| Plant 4 | Upper Root | Positive | 26.78 | 23.17 | 26.05 | 4316.22 |
| Plant 1 | Upper Stem Clone | Positive | 25.54 | 15.66 | 17.04 | 890522.17 |
| Plant 1 | Upper Stem Clone | Positive | 25.39 | 15.12 | 17.13 | 1306380.74 |
| Plant 1 | Upper Stem Clone | Low Positive | 26.26 | 31.26 | 34.48 | 13.86 |
| Plant 1 | Upper Stem Clone | Positive | 25.12 | 16.85 | 18.09 | 382730.67 |
| Plant 1 | Upper Stem Clone | Negative | 26.69 | NaN | NaN | Not Detected |
| Plant 1 | Upper Stem Clone | Low Positive | 26.32 | 35.07 | NaN | 0.93 |
| Plant 2 | Upper Stem Clone | Negative | 25.51 | NaN | NaN | Not Detected |
| Plant 2 | Upper Stem Clone | Negative | 24.09 | NaN | NaN | Not Detected |
| Plant 2 | Upper Stem Clone | Negative | 26.48 | NaN | NaN | Not Detected |
| Plant 2 | Upper Stem Clone | Negative | 26.14 | NaN | NaN | Not Detected |
| Plant 2 | Upper Stem Clone | Negative | 25.23 | NaN | NaN | Not Detected |
| Plant 2 | Upper Stem Clone | Low Positive | 25.96 | NaN | NaN | Not Detected |
| Plant 3 | Upper Stem Clone | Negative | 25.49 | NaN | NaN | Not Detected |
| Plant 3 | Upper Stem Clone | Low Positive | 26.31 | 33.16 | NaN | 3.60 |
| Plant 3 | Upper Stem Clone | Low Positive | 25.04 | 31.66 | NaN | 10.44 |
| Plant 3 | Upper Stem Clone | Low Positive | 25.55 | 30.98 | 43.89 | 16.91 |
| Plant 3 | Upper Stem Clone | Low Positive | 24.46 | 33.5 | NaN | 2.83 |
| Plant 3 | Upper Stem Clone | Low Positive | 25.94 | 31.16 | NaN | 14.88 |

Table 3. (below) Amplification results among samples taken from plants with no known infection by HLVd

| Plant | Tissue | Result | Control CT Value | HLVd Target 1 CT Value | HLVd Target 2 CT Value | Viroid Load (copies/ μ l) |
|----------------------|-------------|--------------|------------------|------------------------|------------------------|-------------------------------|
| Plant 5 | Root | Negative | 26.03 | NaN | NaN | Not Detected |
| Plant 5 | Old Petiole | Negative | 27.16 | NaN | NaN | Not Detected |
| Plant 5 | Old Petiole | Negative | 26.50 | NaN | NaN | Not Detected |
| Plant 5 | Old Petiole | Negative | 27.34 | NaN | NaN | Not Detected |
| Plant 5 | Old Petiole | Negative | 27.23 | NaN | NaN | Not Detected |
| Plant 5 | New Petiole | Negative | 28.22 | NaN | NaN | Not Detected |
| Plant 5 | New Petiole | Negative | 28.41 | NaN | NaN | Not Detected |
| Plant 5 | New Petiole | Negative | 27.86 | NaN | NaN | Not Detected |
| Plant 5 | New Petiole | Negative | 30.06 | NaN | NaN | Not Detected |
| Plant 5 | New Petiole | Negative | 28.31 | NaN | NaN | Not Detected |
| Plant 6 | Root | Negative | 25.23 | NaN | NaN | Not Detected |
| Plant 6 | Old Petiole | Negative | 31.0 | NaN | NaN | Not Detected |
| Plant 6 | Old Petiole | Negative | 28.86 | NaN | NaN | Not Detected |
| Plant 6 | Old Petiole | Negative | 28.98 | NaN | NaN | Not Detected |
| Plant 6 | Old Petiole | Negative | 29.43 | NaN | NaN | Not Detected |
| Plant 6 | New Petiole | Negative | 32.27 | NaN | NaN | Not Detected |
| Plant 6 | New Petiole | Negative | 30.44 | NaN | NaN | Not Detected |
| Plant 6 | New Petiole | Negative | 32.88 | NaN | NaN | Not Detected |
| Plant 6 | New Petiole | Negative | 33.28 | NaN | NaN | Not Detected |
| Plant 6 | New Petiole | Negative | 33.30 | NaN | NaN | Not Detected |
| Control | - | Not Infected | 27.11 | NaN | NaN | Not Detected |
| Control | - | Infected | 28.10 | 15.47 | 18.29 | 1019064.54 |
| Non-Template Control | - | NTC | NaN | NaN | NaN | Not Detected |

Table 4. (below) Amplification results from plants assayed in the follow-up study. Only plants with at least one sample testing positive for HLVd are shown.

| Plant | Tissue | Result | Control CT Value | HLVd Target 1 CT Value | HLVd Target 2 CT Value | Viroid Load (copies/ μ l) |
|---------|--------------|----------|------------------|------------------------|------------------------|-------------------------------|
| Plant 1 | Petiole | Positive | 25.05 | 17.22 | 17.69 | 294348.14 |
| Plant 1 | Petiole | Negative | 26.93 | NaN | NaN | NaN |
| Plant 1 | Petiole | Positive | 22.31 | 18.48 | 19.22 | 120374.98 |
| Plant 1 | Petiole Pool | Positive | 26.73 | 21.04 | 22.33 | 19568.53 |
| Plant 1 | Root | Positive | 26.69 | 22.44 | 24.05 | 7245.79 |
| Plant 1 | Root | Positive | 25.09 | 24.38 | 25.58 | 1828.89 |
| Plant 1 | Root | Positive | 27.13 | 23.95 | 25.08 | 2481.49 |
| Plant 1 | Root Pool | Positive | 25.05 | 17.22 | 17.69 | 294348.14 |
| Plant 2 | Petiole | Negative | 27.25 | NaN | NaN | NaN |
| Plant 2 | Petiole | Negative | 25.92 | NaN | NaN | NaN |
| Plant 2 | Petiole | Negative | 26.43 | NaN | NaN | NaN |
| Plant 2 | Petiole Pool | Negative | 24.39 | NaN | NaN | NaN |
| Plant 2 | Root | Positive | 25.08 | 23.60 | 25.10 | 3181.12 |
| Plant 2 | Root | Positive | 28.19 | 23.32 | 24.18 | 3880.38 |
| Plant 2 | Root | Positive | 27.57 | 26.23 | 30.03 | 492.07 |
| Plant 2 | Root Pool | Positive | 22.89 | 20.96 | 23.68 | 20711.61 |
| Plant 3 | Petiole | Positive | 25.07 | 17.74 | 18.18 | 203516.75 |
| Plant 3 | Petiole | Positive | 23.16 | 16.48 | 17.18 | 497651.41 |
| Plant 3 | Petiole | Negative | 26.01 | NaN | NaN | NaN |
| Plant 3 | Petiole Pool | Positive | 23.08 | 15.83 | 16.53 | 789317.39 |
| Plant 3 | Root | Positive | 28.73 | 21.58 | 22.85 | 13339.31 |
| Plant 3 | Root | Positive | 30.02 | 24.20 | 26.44 | 2078.09 |
| Plant 3 | Root | Positive | 29.48 | 23.36 | 24.98 | 3771.79 |
| Plant 3 | Root Pool | Positive | 26.49 | 24.03 | 25.72 | 2344.53 |

Contact the TUMI Team

With a background in human diagnostics, the scientists at TUMI Genomics thoroughly validate their procedures and believe in the power of collaboration in the scientific community.

Contact TUMI Genomics

Website: www.tumigenomics.com

Email: sales@tumigenomics.com

Phone: (720) 807-8864

Connect with Us online:

Instagram: @TumiGenomics

LinkedIn: TUMI Genomics

Twitter: @TumiGenomics